

ADULT & PEDIATRIC DISPOSABLE BITE BLOCKS



Product Description

The unique ergonomic design of the Clinical Choice Disposable Bite Block provides comfort for the patient as well as ease of oral access for the assistant. Our Adult Bite Block is 60fr and our Pediatric Bite Block is 42fr. Each bite block is latex free, has a dental rim, and the latex free universal strap comfortably fits each patient.

- Packaged individually
- 50/bag
- Non-sterile

Features

- Ample port sizes for scope and dilator. Reduces endoscope friction, allows suction.
- Latex Free safe for use on all patients.
- Universal Adjustable Strap comfortably fits the needs of all patients.
- Dental Rim maintains position, reduces gum damage.
- Solid Construction protects both the patient and valuable equipment.
- Single Patient Use reduces risk of cross contamination.

CCBB-L



- Adult/Large Bite Block
- 20x27 mm - 60 Fr
- Latex Free
- Adjustable Strap
- Call for pricing - 50/bag

CCBB-S



- Pediatric/Small Bite Block
- 14x20 mm - 42 Fr
- Latex Free
- Adjustable Strap
- Call for pricing - 50/bag

Vendor Set Up Info

Federal ID No.
56-2245673

Website
www.clinicalchoice.com

Phone
877.572.9985

Fax
877.572.9986

Mailing Address
5603 W. Friendly Ave. Ste B
PMB 255
Greensboro, NC 27410

Remit to Address
PO Box 890368
Charlotte, NC 28289-0368

Place Your Order

Orders may be placed by:

Phone
877.572.9985

Fax
877.572.9986

Email
orders@clinicalchoice.com



Contact Customer Service with any questions you may have.