

## Simple, Pre-measured Pocket Saves You Time and Money While Improving the Cleaning Process

OneTAB by ClinicalChoice is a family of concentrated powder detergents designed to improve the manual cleaning of endoscopes and surgical instruments.



### Special Formula

Each packet of OneTAB provides a high-dose concentration of three enzymes and detergents that delivers maximum cleaning results during those crucial first minutes of the manual cleaning process. Our proprietary formula ensures performance against a broad spectrum of cleaning challenges.

### Unique Packaging

OneTAB enzymes and detergents are packaged in air-tight packets in their purest form to preserve freshness until ready for use. Each pre-measured packet guarantees correct concentrations of enzymes and detergents every time. This eliminates the need for exact ratios and water volumes required by bulky liquid detergents.

### Better Outcomes. Better Compliance.

The science behind OneTAB delivers the ultimate in convenience, flexibility and streamlined workflow helping you increase endoscope and instrument turnovers. The OneTAB product family offers applications for all of your manual cleaning needs. OneTAB is designed to break down, lift and remove:

- w Blood
- w Protein
- w Carbohydrates
- w Polysaccharides
- w Oils
- w Fats
- w Synthetic Lipids
- w Biofilm
- w Simethicone
- w Cellulose-based adhesives

# OneTAB® ADDITIONAL KEY BENEFITS

			
<h3>COMPLIANT</h3> <ul style="list-style-type: none"> <li>❖ One packet per sink, basin, or ultrasonic cleaner</li> <li>❖ No temperature requirement</li> <li>❖ Neutral pH, low foaming</li> <li>❖ Eliminates test strips</li> <li>❖ Works in hard, soft &amp; deionized waters</li> <li>❖ Safe for scopes &amp; instruments</li> <li>❖ Better storage &amp; transport temperature thresholds                     <ul style="list-style-type: none"> <li>▶ Heat up to 150° F</li> <li>▶ Cold - N/A</li> </ul> </li> </ul>	<h3>EFFICIENT</h3> <ul style="list-style-type: none"> <li>❖ Peak enzyme and detergent efficiency</li> <li>❖ Cleans porous and non-porous surfaces</li> <li>❖ Detergent-only formula available for sensitive and ophthalmic instruments</li> </ul>	<h3>SAVINGS</h3> <ul style="list-style-type: none"> <li>❖ Costs less than bottled detergents</li> <li>❖ Reduces shipping costs</li> <li>❖ Saves valuable storage space</li> <li>❖ One container weighing 5 lbs. replaces (8) 1-gallon bottles of liquid detergent weighing 64 lbs.</li> </ul>	<h3>GREEN</h3> <ul style="list-style-type: none"> <li>❖ 100% biodegradable</li> <li>❖ Allergy friendly</li> <li>❖ Water soluble</li> <li>❖ No dyes or fragrances that can prohibit enzyme efficacy</li> </ul>

## OneTAB® Family of Powder Detergent Products

OneTAB Product	Description	Application	Quantity	Order Number
<b>SINKSIDE™ 3.8-26.5 Litre</b>	Triple Enzymatic Detergent for Manual Cleaning of Scopes and Instruments	Standard sink capacity	200 packets/CNT	OTSS-P
<b>SINKSIDE™ 26.5-53 Litre</b>	Triple Enzymatic Detergent for Manual Cleaning of Scopes and Instruments	Large sink capacity	100 packets /CNT	OTSS-LS
<b>FLEX™ 500 ML- 3.8 Litre</b>	Non-Enzymatic Triple Detergent for Manual Cleaning of Scopes and Instruments. Safe for Ophthalmic Instruments	Bedside, pre-si sink and ultrasonic cleaning.	200 packets/CNT	OT-FLEX
<b>SUPER ULTRASONIC™ 3.8-26.5 Litre</b>	Triple Enzymatic Detergent for Ultrasonic Cleaning of Instruments	Ultrasonic cleaning	100 packets/CNT	OTSU-P

Contact us today to request  
**TRIAL Samples of OneTAB®**

@ [freetrial@clinicalchoice.com](mailto:freetrial@clinicalchoice.com) 1-877-572-9985

0 7300 Cessna Drive Greensboro, NC 27409  
5 1-877-572-9985  
0 info@clinicalchoice.com

